PANSERIN™ 413

PANSERIN™ 413 is a medium ready for use for the cultivation of lymphocytes (from full blood). There is a multitude of questions with the subject matter of short-term cultivation of human cells. As human full blood is relatively easily accessible and because of the possibility of preparing peripheral lymphocyte cultures for a short time, quite a lot of methods have been developed in order to keep blood cells and especially lymphocytes in culture. In principle blood cells in the culture die relatively soon, only lymphocytes can be kept in culture over several cell divisions. However, in order to reach a division of the cells which are normally not proliferating, the cells must be stimulated with specific mitogens. These mitogens are mostly plant lectines (phytohaemagglutinin, PHA).

Composition:

Based on Iscove's MEM, trace elements, albumin, cholesterol, soya-lipids and vitamins were added to the medium. A growth factor mixture is also supplied which has to be added to the medium shortly before the cultivation.

Application:**1. Isolation of lymphocytes from full blood by means of density gradient centrifugation**

- Mix heparinised blood 1:1 with PBS and add it into a centrifuge tube which has been filled with lymphocyte separating medium (PancoL density 1,077 g/ml) before: pipette carefully in order to avoid a phase mixture!
- Centrifuge the gradient at 400 x g for 30 minutes at ambient temperature (brake of the centrifuge set on "off"); 4 phases are created
 - top phase plasma
 - opaque whitish bands (lymphocytes)
 - separating medium
 - pellet with erythrocytes and granulocytes
- Suck off the plasma with a pipette and transfer the lymphocytes with another pipette into a new centrifuge glass.
- Wash the lymphocytes with PBS (without Ca, Mg) and centrifuge them at 100 x g for 10 minutes. Repeat washing step.

2. Cultivation and stimulation of the lymphocytes

- Resuspend lymphocytes in PANSERIN™ 413.
- For the stimulation of the lymphocyte proliferation adjust the cells to a cell density of approx. 1×10^5 and add phytohaemagglutinin in a concentration of 1 - 5 mg/ml; the incubation time is 48 to 72 hours, depending on the kind and origin of the lymphocytes and depending on the further use.
- Further cultivation in PANSERIN™ 413. After approx. 14 days a restimulation is necessary.

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BIOTECH GmbH