
PANSERIN™ PX10


PANSERIN™ PX10 is a serumfree medium ready for use for the production of monoclonal antibodies. Due to the optimised composition the cells are stimulated to proliferate even with low seed densities.

SP2/O-Ag-14 in PANSERIN™ PX10
SP2/O-Ag-14 in PANSERIN™ PX10

SP2/O-Ag-14 have been transferred from a culture containing serum (10 % FCS in RPMI 1640) directly into PANSERIN™ PX10 and cultivated in the incubator at 37 °C and with 5 % CO₂-fumigation for 5 days. Seed 1000 cells/ml. The cells grow very well without a long adaption phase, comparable to a culture containing serum (picture on the right).


PANSERIN™ PX10 distinguishes itself, beside the excellent growth-stimulating properties, by its cloning capacity. Conventional serumfree systems often need lengthy and timeconsuming adaption steps and cell seeds up to 10⁵ cells/ml. By contrast, most clones can be directly transferred into PANSERIN™ PX10. With PANSERIN™ PX10, clones can even be produced without great difficulties.

Cloning of SP2/O-Ag-14 in PANSERIN™ PX10
Klonbildung von SP2/O-Ag-14 in PANSERIN™ PX10

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
Cloning of SP2/O-Ag-14 in RPMI 1640 with 20 % FCS
Klonbildung v. SP2/O-Ag-14 in RPMI 1640 m. 20 % FCS

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Am Gewerbepark 13, D-94501 Aidenbach
 Phone +49 (0) 85 43 - 60 16-30, Fax +49 (0) 85 43 - 60 16-49
 e-mail: info@pan-biotech.de, http://www.pan-biotech.de


PAN™
 BIOTECH GmbH

Composition:

PANSERIN™ PX10 contains cleaned proteins, lipids, salts, amino acids, trace elements and hormones in an optimised prescription. It doesn't contain any undefined hydrolysates or lysates (e.g. peptones). L-glutamine is already included in the medium. Antibiotic is not included in the formulation.

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PAN[™]
BIOTECH GmbH

Monoclonal Antibodies

Antibodies are proteins specifically binding to structures res. substances which are called antigens. They are part of the defense mechanisms of the immune system and are produced by the leucocytes (B-lymphocytes) as a reaction on the respective specific antigen.

The human humoral immune system is able to produce a tremendous number of antibodies (approx. 10^8) with different antigen binding specificity. They serve to identify foreign body antigens (bacteria, virus).

This diversity is brought about by a relatively small gene repertoire and a recombination of V-, J- and D-sequences takes effect on the variety of antibodies of the variable regions.

The work of Köhler and Milstein set a landmark for the production of monoclonal antibodies. The conventional way of production used to bring about only polyclonal antibodies in small quantities but the specificity and presentable quantity for experimental work were made possible only by the hybridisation of cells (not indefinitely divisible, e.g. spleen cells) from immunized animals and the indefinitely divisible myeloma cells. Such hybrid cells can be multiplied further on in the culture and then they produce a clone (offspring of identical cells) by which monoclonal antibodies are produced.

The field of application and the distribution of monoclonal antibodies reach from a therapeutical application (fight against cancer) to diagnostic applications (ELISA) and are still increasing.

Application:

In many cases a serumfree cultivation can be carried out without time-consuming adaption steps. With critical clones, reduce the serum slowly (see instruction PANSERIN™ 401).

- Warm up PANSERIN™ PX10 to 37 °C.
- Transfer hybridoma directly into PANSERIN™ PX10. In most cases a cell number of 1000 cells/ml is sufficient.
- Incubate the cells in the usual way in the CO₂-incubator at 37 °C (5 % CO₂-fumigation) res. multiply them in the bioreactor.
- Extract monoclonal antibodies from the supernatant.

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
PAN[™]
BIOTECH GmbH


Typical growth curve of SP2/0-Ag-14 in PANSERIN™ PX10

Sp2/0-Ag-14 cells have been transferred from a culture containing serum (RPMI 1640 with 10 % FCS) directly into PANSERIN™ PX10. Seed 1000 Zellen/ml. As a comparison Sp2/0-Ag-14 in RPMI 1640 with 10 % FCS.

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SP2/0-Ag-14 cells were used from a culture containing serum directly in a cloning efficiency test with PANSERIN™ PX10

The absolute cloning efficiency was 75 % (77 % with 20 % FCS) the relative cloning efficiency 97 % (100 % with 20 % FCS)

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