

NEW

PANSERIN H4000

proteinfree
cell culture medium
for hybridoma cells

ready to use
without any supplements



PANSERIN H 4000

for Hybridoma
proteinfree
with: L-Glutamine
sterile filtered

Cat.No.: P04-714000
Lot.No.: 200707
Volume: 500 ml
Store at: +4 °C
Expiry: 07/2008

for cell culture

PANTM
BIOTECH GmbH

**PAN
BIOTECH**
Just ask us!




www.pan-biotech.de
**THE BEST
IN
CELL CULTURE**

See the back of the
page for an overview
of all the advantages

PANTM
BIOTECH GmbH



PANSERIN H4000



ADVANTAGES:

- free from human and animal components**
- short adaption phase** (possibly no adaption required)
 - suitable for a wide range of hybridomas
- excellent growth properties**
- very high antibody production**
- simple and **rapid cleaning of antibodies**
- also suitable for spinner cultures**

ADVANTAGES of a serumfree cell culture:

- lower risk** with regard to contamination
 - with bacteria, fungi or virus
- better defined and reproducible formulations**
 - allow more convincing and comparable research results
- time-consuming **batch tests are not necessary**
- elimination of a source for possible infectious agents**
- facilitation of the cleaning** of the end products
 - Fulfilment of legal conditions
- for the production of medical products**
- reduction of contaminations**
 - of the end products by culture residues

**PAN
BIOTECH**
Just ask us!

Am Gewerbepark 13
D-94501 Aidenbach
Tel. +49(0)8543/60 16-30
Fax +49 (0)8543/60 16-49
info@pan-biotech.de

www.pan-biotech.de
**THE BEST
IN
CELL CULTURE**

PANTM
BIOTECH GmbH

